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# Summary of Disk Diffusion Susceptibility Testing Conditions


Organism Groups	Test Medium (Agar)	Inoculum Size(CFU/mL)	Incubation Conditions	Incubation Duration
Enterobacterales	Mueller-Hinton Agar (MHA)	Swab from $1.5 \times 10^8$ suspension	35°C; room air	16–18 h
<i>Pseudomonas aeruginosa</i>	MHA	Swab from $1.5 \times 10^8$ suspension	35°C; room air	16–18 h
Enterococci	MHA	Swab from $1.5 \times 10^8$ suspension	35°C; room air	16–18 h (24 h for vancomycin R <sup>a</sup> )
Staphylococci (to detect methicillin-resistance)	MHA (add 2% NaCl for meth-R)	Swab from $1.5 \times 10^8$ suspension	30°–35°C; air (temp >35°C may mask meth-R)	16–18 h (24 h for methicillin and vancomycin R <sup>a</sup> )
<i>Streptococcus pneumoniae</i> and other streptococci	MHA plus 5% sheep blood	Swab from $1.5 \times 10^8$ suspension	35°C; 5%–7% carbon dioxide (CO <sub>2</sub> )	20–24 h
<i>Haemophilus influenzae</i>	<i>Haemophilus</i> test medium	Swab from $1.5 \times 10^8$ suspension	35°C; 5%–7% CO <sub>2</sub>	16–18 h
<i>Neisseria gonorrhoeae</i>	GC agar plus supplements	Swab from $1.5 \times 10^8$ suspension	35°C; 5%–7% CO <sub>2</sub>	20–24 h

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
**FIGURE 13-1** ■ The tube on the left is a McFarland 0.5 turbidity standard. The tube on the right is a test bacterial suspension that has a turbidity greater than that of the McFarland standard; it is more difficult to see the black lines through the test suspension than through the McFarland standard. Sterile saline or broth must be added to the test suspension to dilute it until the turbidity matches that of the McFarland standard.

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within **15 minutes** after adjusting the turbidity, a **sterile cotton swab** is dipped into the suspension.

The swab should be rotated several times and pressed firmly on the inside wall of the tube above the fluid level to remove excess inoculum from the swab.

The Mueller-Hinton agar plate is inoculated by streaking the swab

This procedure is repeated by streaking two more times, rotating the plate approximately 60° each time to ensure an even distribution of inoculum.

As a final step, the rim of the agar is swabbed

The lid may be left ajar for 3 to 5 minutes, but no more than 15 minutes, to allow for any excess surface moisture to be absorbed before applying the drug-impregnated disks.

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**GROUP A  
PRIMARY TEST  
AND REPORT**

**GROUP B<sup>e</sup>  
PRIMARY TEST  
REPORT SELECTIVELY**

**GROUP C<sup>f</sup>  
SUPPLEMENTAL  
REPORT SELECTIVELY**

**GROUP U  
SUPPLEMENTAL  
FOR URINE ONLY**

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## Anatomic Distribution of Some Common Antibacterial Agents

	Serum or Blood <sup>a</sup>	Cerebrospinal Fluid	Urine
<b>Ampicillin</b>	+	+	+
<b>Ceftriaxone</b>	+	+	+
<b>Meropenem</b>	+	+	+
<b>Vancomycin</b>	+	±	+
<b>Ciprofloxacin</b>	+	±	+
<b>Gentamicin</b>	+	-	+
<b>Clindamycin</b>	+	-	-
<b>Norfloxacin</b>	-	-	+
<b>Nitrofurantoin</b>	-	-	+

<sup>a</sup> Serum or blood represents a general anatomic distribution.

+, Therapeutic levels generally achievable at that site;

±, therapeutic achievable levels moderate to poor;

-, therapeutic levels generally not achievable at that site.

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

FIG. 11.7 Examination of a disk diffusion plate by transmitted and reflected light.

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
**FIGURE 13-10** ■ Routine disk diffusion tests are examined by placing the plate on or 2 to 3 inches above a black, nonreflecting surface. Reflected light is used to illuminate the plate.

## Modify methods for fastidious bacteria


If blood was added to the agar base (as with streptococci), the zones are measured from the upper surface of the agar illuminated with reflected light, with the cover removed.

**FIGURE 13-11** ■ Disk diffusion tests for staphylococci with oxacillin (or methicillin or nafcillin) and vancomycin and for enterococci with vancomycin are examined by holding the plate up to a light source (transmitted light) for zone examination. Any growth within the zone is significant.

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

**Figure 12-9** Individual bacterial colonies within a more obvious zone of inhibition (*arrows*). This could indicate inoculation with a mixed culture. However, emergence of resistant mutants of the test isolate is a more likely reason for this growth pattern.


**FIGURE 13-13** ■ The oxacillin zone for heteroresistant oxacillin-resistant *Staphylococcus aureus* often shows a haze of growth within the zone of inhibition. This haze is significant, and the isolate here is oxacillin resistant.


**FIGURE 13-9** ■ *Escherichia coli* tested by the disk diffusion method. The lawn of growth following overnight incubation shows individual colonies, representing unsatisfactory growth. The most likely explanation for the scanty growth is the use of an inoculum that either is too light or contains too many nonviable cells, resulting in larger than normal zones and potentially false susceptible results.

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## Appendix E. Dosage Regimens Used to Establish Susceptible or Susceptible-Dose Dependent Breakpoints

The evolving science of pharmacokinetics-pharmacodynamics has become increasingly important in recent years in determining minimal inhibitory concentration (MIC) breakpoints. Recently approved susceptible or susceptible-dose dependent (SDD) breakpoints for a number of agents have been based on a specific dosage regimen(s); these dosage regimens are listed in the table below. Proper application of the breakpoints necessitates drug exposure at the site of infection that corresponds to or exceeds the expected systemic drug exposure at the dose listed in adult patients with normal renal function. This information should be shared with pharmacists, infectious diseases staff, and others making dosing recommendations for the institution.

Antimicrobial Agent	Breakpoints and Interpretive Categories			
	Susceptible		SDD	
	MIC	Dose	MIC	Dose
<b>Table 2A. Enterobacterales</b>				
Azithromycin ( <i>Salmonella enterica</i> ser. Typhi)	≤16 µg/mL	500 mg administered daily	N/A	
Aztreonam	≤4 µg/mL	1 g administered every 8 h	N/A	
Cefazolin	≤2 µg/mL	2 g administered every 8 h	N/A	
Ceftaroline	≤0.5 µg/mL	600 mg administered every 12 h	N/A	
Cefepime	≤2 µg/mL	1 g administered every 12 h	4 µg/mL	1 g administered every 8 h or 2 g administered every 12 h
			8 µg/mL or zone diameter: 19–24 mm	2 g administered every 8 h  (Because it is not possible to correlate specific zone diameters with specific MICs, an isolate with a zone diameter in the SDD range should be treated as if it might be an MIC of 8 µg/mL.)

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Example of Breakpoints and Interpretive Categories as Used in Table 2

Antimicrobial Agent	Disk Content	Interpretive Categories and Zone Diameter Breakpoints, nearest whole mm			Interpretive Categories and MIC Breakpoints, $\mu\text{g}/\text{mL}$		
		S	I <sup>*</sup>	R	S	I <sup>*</sup>	R
X	30 $\mu\text{g}$	$\geq 20$	15–19	$\leq 14$	$\leq 4$	8–16	$\geq 32$
Y	–	–	–	–	$\leq 1$	2	$\geq 4$
Z	10 $\mu\text{g}$	$\geq 16$	–	–	$\leq 1$	–	–

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Organism	Resistance Phenotype Detected
Streptococcus pneumoniae	penicillin – R Vancomycin – NS
Streptococcus, $\beta$ -hemolytic group	penicillin – NS Vancomycin – NS
Streptococcus, viridans group	Vancomycin – NS
Enterococcus spp.	Linezolid ,Vancomycin, High-level aminoglycoside – R
Enterococcus faecalis	Linezolid ,penicillin, High-level aminoglycoside – R
Enterococcus faecium	Linezolid, High-level aminoglycoside – R
Any Organism	Resistance to all agents routinely tested

Organism	Resistance Phenotype Detected
Acinetobacter baumannii	Carbapenem – I or R
Pseudomonas aeruginosa	Colistin – I or R    Carbapenem – I or R cuncurrent Gentamicin,Tobramycin and Amikacin - R
Stenotrophomonas maltophilia	Cotrimoxazole - I or R
Any Organism	Resistance to all agents routinely tested

Organism	Resistance Phenotype Detected
Staph.aureus	Linezolid,Cefoxitin - R      Vancomycin - R
Staph.coagulase neg	Linezolid - R                      Vancomycin - R
Any Organism	Resistance to all agents routinely tested

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Screen Test	Vancomycin MIC $\geq 8 \mu\text{g/mL}$	
Test method	Agar Dilution	Agar Dilution
Organism group	<i>S. aureus</i>	<i>Enterococcus</i> spp.
Medium	BHI agar	BHI <sup>a</sup> agar
Antimicrobial concentration	6 $\mu\text{g/mL}$ vancomycin	6 $\mu\text{g/mL}$ vancomycin
Inoculum	Colony suspension to obtain 0.5 McFarland turbidity  Preferably, using a micropipette, spot a 10- $\mu\text{L}$ drop onto agar surface. Alternatively, using a swab dipped in the suspension and the excess liquid expressed, spot an area 10–15 mm in diameter or streak a portion of the plate.	1–10 $\mu\text{L}$ of a 0.5 McFarland suspension spotted onto agar surface. Alternatively, using a swab dipped in the suspension and the excess liquid expressed, spot an area 10–15 mm in diameter or streak a portion of the plate.
Incubation conditions	35°C $\pm$ 2°C; ambient air	35°C $\pm$ 2°C; ambient air
Incubation length	24 hours	24 hours
Results	Examine carefully with transmitted light for >1 colony or light film of growth.  >1 colony = Presumptive reduced susceptibility to vancomycin	>1 colony = Presumptive vancomycin resistance
Additional testing and reporting	Perform a vancomycin MIC using a validated MIC method to determine vancomycin MICs on <i>S. aureus</i> that grow on BHI–vancomycin screening agar.  Testing on BHI–vancomycin screening agar does not reliably detect all vancomycin-intermediate <i>S. aureus</i> strains. Some strains for which the vancomycin MICs are 4 $\mu\text{g/mL}$ will fail to grow.	Perform vancomycin MIC on <i>Enterococcus</i> spp. that grow on BHI–vancomycin screening agar and test for motility and pigment production to distinguish species with acquired resistance (eg, <i>vanA</i> and <i>vanB</i> ) from those with intrinsic, intermediate-level resistance to vancomycin (eg, <i>vanC</i> ), such as <i>Enterococcus gallinarum</i> and <i>Enterococcus casseliflavus</i> , which often grow on the vancomycin screen plate. In contrast to other enterococci, <i>E. casseliflavus</i> and <i>E. gallinarum</i> with vancomycin MICs of 8–16 $\mu\text{g/mL}$ (intermediate) differ from vancomycin-resistant enterococci for infection <b>prevention</b> purposes.

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## Screening Tests

Screening Test	Organisms	Test Description	When to Perform Confirmatory Test	Confirmatory Test
Vancomycin agar screen	<ul style="list-style-type: none"> <li><i>S. aureus</i></li> <li><i>Enterococcus</i> spp.</li> </ul>	Agar dilution; BHI with 6 µg/mL vancomycin	If screen positive	Vancomycin MIC

Table 2D. *Enterococcus* spp. (Continued)

Test/Report Group	Antimicrobial Agent	Disk Content	Interpretive Categories and Zone Diameter Breakpoints, nearest whole mm			Interpretive Categories and MIC Breakpoints, µg/mL				Comments
			S	I	R	S	SDD	I	R	
<b>GLYCOPEPTIDES</b>										
B	Vancomycin	30 µg	≥17	15–16	≤14	≤4	–	8–16	≥32	<p>(10) When testing vancomycin against enterococci, plates should be held a full 24 hours for accurate detection of resistance. Zones should be examined using transmitted light; the presence of a haze or any growth within the zone of inhibition indicates resistance. Organisms with intermediate zones should be tested by an MIC method as described in M07.<sup>3</sup> For isolates for which the vancomycin MICs are 8–16 µg/mL, perform biochemical tests for identification as listed under the “Vancomycin MIC ≥ 8 µg/mL” test found in Table 3G.</p> <p>See general comment (4) and comment (8).</p>

**GLYCOPEPTIDES**

(20) MIC tests should be performed to determine the susceptibility of all isolates of staphylococci to vancomycin. The disk test does not differentiate vancomycin-susceptible isolates of *S. aureus* from vancomycin-intermediate isolates, nor does the test differentiate among vancomycin-susceptible, -intermediate, and -resistant isolates of *Staphylococcus* spp. other than *S. aureus*, all of which give similar size zones of inhibition.

B	Vancomycin	<i>S. aureus</i>	-	-	-	-	-	≤2	-	4-8	≥16	<p>(21) For <i>S. aureus</i>, vancomycin-susceptible isolates may become vancomycin intermediate during the course of prolonged therapy.</p> <p>(22) Send any <i>S. aureus</i> for which the vancomycin is ≥8 µg/mL to a referral laboratory. See Appendix A.</p> <p>Also refer to Table 3F for <i>S. aureus</i>, Subchapter 3.12 in M07,<sup>3</sup> and Subchapter 3.9 in M02.<sup>1</sup></p>
		<i>Staphylococcus</i> spp. other than <i>S. aureus</i>	-	-	-	-	-	≤4	-	8-16	≥32	<p>See comment (19).</p> <p>(23) Send any <i>Staphylococcus</i> spp. other than <i>S. aureus</i> for which the vancomycin MIC is ≥32 µg/mL to a referral laboratory. See Appendix A.</p> <p>See also Subchapter 3.12 in M07<sup>3</sup> and Subchapter 3.9 in M02.<sup>1</sup></p>

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

❖ for S.aureus and S.lugdunensis

$\geq 22$  mm.... (S) and  $\leq 21$  mm .... (R)


❖ for coagulase-negative staphylococci :

$\geq 25$  mm .... (S) and  $\leq 24$  mm .... (R)

## Supplemental Methods for Detection of Antimicrobial Resistance

Test	Purpose (Bacteria/Antimicrobial-R)	Conditions
<b>Cefoxitin to detect <i>mecA</i>-mediated oxacillin resistance</b>	Staphylococci ( <i>S. aureus</i> and <i>S. lugdunensis</i> ) and coagulase-negative staphylococci (CoNS)/oxacillin-resistance	Disk diffusion: MHA and 30 µg cefoxitin disk Incubation: 33°C–35°C, 16–18 h: <i>S. aureus</i> and <i>S. lugdunensis</i> ; 24 h: CoNS


**FIGURE 13-14** ■ D-zone testing of *Staphylococcus aureus* that has erm-mediated inducible clindamycin resistance. The positive reaction is noted by a flattening of the zone around the clindamycin disk in the area where there has been diffusion of both erythromycin and clindamycin molecules.

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## Supplemental Methods for Detection of Antimicrobial Resistance

Test	Purpose (Bacteria/Antimicrobial-R)	Conditions
D-zone test	<p><i>S. aureus</i>, <i>S. lugdunensis</i>, CoNS, <i>Streptococcus pneumoniae</i>, and beta-hemolytic streptococci (e.g., group B)/inducible clindamycin resistance</p> <p>Differentiation of inducible iMLS<sub>B</sub> resistance (by <i>ermA</i> or <i>ermC</i> genes) versus erythromycin resistance by efflux mechanism (i.e., <i>msrA</i> gene)</p>	<p>Medium: 5% sheep blood agar or MHA (for staphylococci) or MHA plus 5% sheep blood (for streptococci)</p> <p>Antimicrobials: 15 µg erythromycin (E) and 2 µg clindamycin (Cd) (for staphylococci, space 15–26 mm apart; for streptococci, space 12 mm apart)</p> <p>Incubation:            staphylococci: 35°C, air, 16–18 h;            streptococci: 35°C, 5% CO<sub>2</sub>, 20–24 h</p>

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**Figure 59-4** Disk induction or “D” test. The *Staphylococcus aureus* isolate on the left demonstrates erythromycin (E) resistance but susceptibility to clindamycin (CC). The *S. aureus* isolate on the right demonstrates inducible clindamycin resistance with blunting of the zone of inhibition around the clindamycin disk adjacent to the erythromycin disk.

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# *Nitrocefin disk test (Chromogenic cephalosporin method)*

## *Penicillin Disk Zone-Edge Test*

**Figure 1. Positive for  $\beta$ -Lactamase Detection**

**Figure 2. Negative for  $\beta$ -Lactamase Detection.**

# *Screening for ESBL producers*

- for screening for ESBL production by
  - klebsiellae
  - Escherichia coli
  - Proteus mirabilis

***For K.pneumonia, K.oxytoca  
and E. coli***


- Cefpodoxime < 17mm
- ceftazidime < 22mm
- aztreonam < 27mm
- Cefotaxime < 27mm
- Ceftriaxone < 25mm

***For Proteus mirabilis***

- Cefpodoxime <22 mm
- Ceftazidime <22mm
- Cefotaxime <27mm

# *ESBL Confirmatory Tests*

Disc with  
cephalosporin +  
clavulanic acid

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


# Modified Hodge Test

*E. coli* ATCC® 25922

Inhibition of *E. coli* ATCC® 25922 by ertapenem

Enhanced growth of *E. coli* ATCC® 25922. Carbapenemase produced by *K. pneumoniae* ATCC® BAA-1705 inactivated ertapenem that diffused into the media. Thus, there is no longer sufficient ertapenem here to inhibit *E. coli* ATCC® 25922 and an

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Comment: **Vancomycin** susceptibility must be tested with MIC but its resistance reported very rare.

#### Resistant

- 1- Tetracycline
- 2- Penicillin
- 3- Erythromycin
- 4- Oxacillin
- 5- Clindamycin
- 6- Azithromycin
- 7- Clarithromycin

#### Also Resistant based on Oxacillin Resistance:

Amoxicillin, Ampicillin, Cephalexin, Cloxacillin, Cefazolin, Cefuroxime, Dicloxacillin, Cefepime, Amoxicillin-clavulanate, Cefotaxime, Ampicillin-sulbactam, Piperacillin-tazobactam, Cefotetan, Ceftizoxime, Ceftriaxone, Ceftaroline, Moxalactam, Imipenem, Meropenem.

#### Intermediate

No intermediate antibiotic was observed

Comment: **Below antibiotics cannot be used for treatment of *Staphylococci* spp:**  
Nalidixic acid, Aztreonam, and Polymyxin B/Colistin.

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# Urine Culture Result

Patient ID: NAME: ZAHRA SADAT SHIRZADI REC. NO: 1123 AGE: 11 Years

Pathogen(S): *Streptococcus pyogenes* (Group A)

Count: Many (Colonies)

## Antimicrobial Susceptibility Test

### Comment

Routine Susceptibility testing of *Streptococcus pyogenes* isolates is not advised according to Clinical Laboratory Standard Institute: CLSI; 2020.

### Comment

Antibiotics that can be used for *Streptococcus pyogenes* treatment:

1-Penicillin 2-Ampicillin 3-Amoxicillin 4-Amoxicillin-clavulanic acid 5-Ceftizoxime  
6-Cefazolin 7-Ampicillin/sulbactam 8-Cephradine 9-Cephalothin 10-Cefotaxime  
11-Ceftriaxone 12-Cefepim 13-Imipenem 14-Ertapenem 15-Meropenem 16-Cefaclor  
17-Cefdinir 18-Cefprozil 19-Ceftibuten 20-Cefuroxime 21-Cefpodoxime 22-Cephapirin

Patient ID: NAME: ZAHRA SADAT SHIRZADI REC.NO: 1123 AGE: 11 Years

Pathogen(S): *Streptococcus agalactiae* (Group B)

Count: >100,000 CFU/mL

### Antimicrobial Susceptibility Test

#### Comment

Routine Susceptibility testing of urine isolates of *Streptococcus agalactiae* is not advised according to Clinical Laboratory Standard Institute: CLSI; 2017.

#### Comment

Antibiotics that can be used for *Streptococcus agalactiae* treatment:

1-Penicillin 2-Ampicillin 3-Amoxicillin 4-Amoxicillin-clavulanic acid 5-Cefazolin  
6-Ampicillin/sulbactam 7-Cefepim 8-Cephadrine 9-Cephalothin 10-Cefotaxime  
11-Ceftriaxone 12-Ceftizoxime 13-Imipenem 14-Ertapenem 15-Meropenem

#### Comment

recomendations for intrapartum prophylaxi for Gruop B streptococci are penicilin or ampicilin. Although cefazolion is recomended for penicilin-allergic woman at low risk for anaphylaxis, but for those woman at high risk for anaphylaxis may recive clindamycin.

# Urine Culture Result

Patient ID: NAME: ZAHRA SADAT SHIRZADI REC.NO: 1123 AGE: 11 Years

Pathogen(S): *Staphylococcus saprophyticus*

Count: >100,000 CFU/mL

## Antimicrobial Susceptibility Test

### Comment

Routine Susceptibility testing of urine isolates of *Staphylococcus saprophyticus* is not advised according to Clinical Laboratory Standard Institute (CLSI; 2020) because infections of *Staphylococcus saprophyticus* respond to antimicrobial agents concentrations achieved in urine that commonly used to treat acute, uncomplicated urinary tract infections (eg, Nitrofurantoin, Trimethoprim ± Sulfamethoxazole, or a Fluoroquinolone).

### Comment

Antibiotics that can be used for *Staphylococcus saprophyticus* treatment:

1-Nitrofurantoin	2-Trimethoprim	3-Trimethoprim-sulfamethoxazole	4-Ciprofloxacin	
5-Levofloxacin	6-Ofloxacin	7-Moxifloxacin	8-Lomefloxacin	9-Norfloxacin
10-Enoxacin	11-Gatifloxacin	12-Grepafloxacin	13-Sparfloxacin	14-Fleroxacin

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# AST - QC Program

- **Daily Testing**


- ✓ Test all applicable QC strains for *20 or 30 consecutive test days*, and document results

- **Weekly Testing**


- ✓ To convert from daily to weekly QC testing, no more than *one out of 20 or three out of 30* zone diameters for each antimicrobial agent/organism combination may be outside the acceptable zone diameter limits stated

Appendix A. (Continued)

**Disk Diffusion Weekly Quality Control Testing Protocol**

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

**FIGURE 16-11** Etest uses the principle of a predefined antibiotic gradient on a plastic strip to generate an MIC value. It is processed like the disk diffusion. A, Individual antibiotic strips are

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**FIG. 10.6** Cations ( $Mg^{2+}$  and  $Ca^{2+}$ ) and aminoglycosides ( $AG^{++}$ ) compete for the negatively charged binding sites on the outer membrane surface of *P. aeruginosa*. Such competition is an example of the effect that environmental factors (e.g., cation concentrations) can have on the antibacterial activity of aminoglycosides.

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
**Record reaction as:**

4+  
3+  
2+  
1+  
negative

Score and record the results as follows:

4+	All organisms are clumped and the supernatant, or suspending fluid is clear
3+	About 75% agglutination, the supernatant fluid is slightly cloudy
2+	About 50% agglutination, the supernatant fluid is moderately cloudy
1+	About 25% agglutination, the supernatant fluid is cloudy
Tr	Trace amount of agglutination present
Negative	No agglutination apparent and suspension remains homogenous.

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